

Yatee Gupta¹, Yachee Gupta¹, Amar Ranjan², Samrat Mukhopadhyay³


¹Fabiosys Innovations Pvt. Ltd. (IIT Delhi), ²IRCH, AIIMS Delhi, ³Dept. Of Textile and Fiber Engineering, IIT Delhi.

BACKGROUND

- ❖ Patients undergoing surgeries in low-resource settings have a high chance of acquiring common antimicrobial-resistant pathogens like MRSA & VRE from the environment via direct and cross contaminations.
- ❖ It has been reported that antimicrobial-resistant ESKAPE bacteria, viruses, and fungi commonly available in the environment get deposited on fabric surfaces during surgeries in low-resource settings.
- ❖ A proper sterilized environment is needed to enhance the surgical outcome by reducing the chances of infections, mortality, and long-term usage of antibiotics to prevent Surgical Site Infections (SSIs).
- ❖ One of the cheapest ways to provide a sterilized environment during surgery is the use of high-performance antimicrobial surgical supplies like surgical drapes, patient wear, bedsheets, and other medical consumables.

METHODS

- ❖ The commonly used hospital grade cotton, poly-cotton, and polypropylene textile substrates are taken and chemically processed with a metal-free proprietary antimicrobial formulation to provide High-performance Antimicrobial functionality to the fabric, while also making the fabric wash-friendly, comfortable, non-cytotoxic, and biocompatible for human skin.



Input fabric



Required: Arrangement of Textile machineries in a sequence.
Control parameters:

- Temperature
- pH
- Fabric specifications
- Chemical concentrations (optimized)
- Additional Electronic Controls.

Processing Technology



Output fabric

“ Patients have a high chance of catching infections during surgeries. ”



Source: AO Foundation: AO surgery reference.

RESULTS

- ❖ The fabric is rigorously tested for its Antibacterial, Antiviral & Antifungal properties in accordance with AATCC 100, ISO 18184 & AATCC 30 Standards respectively. It is found that >99.9% of the bacteria, viruses, and fungi got destroyed within 30 minutes of contact with fabric.

Test Organisms used		Described by the customer - Fabric: Fabiosys Antimicrobial fabric	
		No. of Swatches Result	Bacterial Reduction (%)
<i>Klebsiella pneumoniae</i> ATCC 4352 - 30 min		5	99.99
<i>Staphylococcus aureus</i> ATCC 6538 - 30 min		5	99.99

Result 1 : The sample showed 99.99 % antibacterial activity against *Staphylococcus aureus* ATCC 6538 and 99.99 % antibacterial activity against *Klebsiella pneumoniae* ATCC 4352 when tested according to AATCC 100-2019 test method.

Results:				
Test Virus: Human Coronavirus HCoV-229E (Surrogate of SARS-CoV-2)				
1: Test Sample: Surface Disinfectant: Fabric: PV: 4000 45 Gam Embossed				
Virus	Replicates	Log of Infectivity titre value Log (Va) (Immediately after inoculation of Control)	Log of Infectivity titre value Log (Vb) (After 30 minutes of contact with Control)	Log of Infectivity titre value Log (Vc) (After 30 minutes of contact with Test Specimen)
Influenza virus suspension: (1.70 × 10 ⁶ PFU / ml)	Set I	5.90	5.85	2.78
	Set II	5.92	5.88	2.85
	Set III	5.92	5.80	2.88
	IgTCD ₅₀ / Average	5.88	5.85	2.84
Reduction value M = Log (Va - Vb)		0.03		-
Log of Antiviral Efficacy Value - Mv		3.01		
Percentage		99.99		
Antiviral Efficacy Value - Mv		99.99		

Test Organism Used		Described by the customer - Fabric: Fabiosys Antimicrobial fabric	
		No. of Swatches Result	Bacterial Reduction (%)
<i>Aspergillus niger</i> ATCC 6275		5	99.99
Zone of inhibition in mm		0 mm	
Growth under fabric		Absent	


Result 1 : The sample showed presence of antifungal activity against *Aspergillus niger* ATCC 6275 when tested according to AATCC 30 Part 3 method.

- ❖ The microbes taken are gram-positive & gram-negative bacteria including MRSA & VRE, enveloped viruses: HCoV-229E and H3N2, and fungi: *Aspergillus Niger*.


Test Organisms used		Described by the customer - Fabric: Fabiosys Antimicrobial fabric	
		No. of Swatches Result	Bacterial Reduction (%)
Methicillin Resistant <i>Staphylococcus aureus</i> -24 hours		5	99.99
Vancomycin Resistant <i>Enterococci</i> - 24 hours		5	99.99

- ❖ The fabric is also tested for its Cytotoxicity via MTT Test method as per ISO 10993-5 which revealed non-toxicity after an exposure of 48 hours, and no skin sensitization or irritation is observed over a duration of 72 hours when tested on New Zealand White Rabbits as per ISO 10993-23 Standard.

Invitro Cytotoxicity- Direct method. Cell line used: L929. Experiment done by: SITRA.



Control



Sample

CONCLUSION

High-performance Antimicrobial surgical supplies can act as an infection control tool against AMR infections.

Use Anti-AMR surgical drapes to prevent SSIs.

Make surgeries safer

LAB TO CLINICS



Surgical Site Infections (SSIs) increase patient mortality, and damage surgeon's, hospital's & medical camp's reputation.

Drape packs are available for every surgery. Contact **partner@fabiosys.com** for a complimentary sample.

Media Coverage:



Acknowledgements

The authors thank Dept. of Biotechnology, Govt. of India to financially support this research & product development. The authors also thank SITRA for independent testing and evaluation of the products.