

Validity and feasibility of a Pan-Lassa rapid diagnostic test for Lassa fever in Abakaliki, Nigeria: a field evaluation



Jelte Elsinga^{1*}, Temmy Sunyoto^{1*}, Letizia di Stefano², Pier Francesco Giorgetti³, Htet Aung Kyi⁴, Chiara Burzio², Ximena Campos Moreno², Chiedozie K Ojide⁵, Nennaya Ajayi⁵, Richard Ewah⁵, Emeka Onwe Ogah⁵, Chioma Dan-Nwafor⁶, Anthony Ahumibe⁶, Chinwe Lucia Ochu⁶, Adebola Olayinka⁶, Sylvie Jonckheere², Pascale Chaillet², Michel van Herp²

¹Luxembourg Operational Research Unit (LuXOR), MSF (Médecins Sans Frontières), Luxembourg, ²MSF, Brussels, Belgium; ³MSF, Abakaliki, Nigeria; ⁴MSF, Abuja, Nigeria; ⁵Alex Ekwueme Federal Teaching Hospital, Abakaliki, Nigeria; ⁶Nigeria Centre for Prevention and Disease Control (NCDC), Abuja, Nigeria

Introduction

- Lassa fever is a viral haemorrhagic fever with few options for diagnosis and treatment;
- Transmitted by rodents – *Mastomys natalensis* and by human (body fluids)
- Endemic in Nigeria, Liberia, Guinea and Sierra Leone
- A point-of-care bedside test diagnosing Lassa fever, adhering to REASSURED criteria, is **not currently available** BUT is urgently needed in west African regions with high Lassa fever burden.
- We aimed to assess the validity and feasibility of a rapid diagnostic test (RDT) to confirm Lassa fever in Nigeria

REASSURED criteria:

Real-time connectivity, ease of specimen collection, affordable, sensitive, specific, user-friendly, rapid and robust, equipment-free or simple, and deliverable to end-users.

Methods

Study design: Prospective study


- Index test:** ReLASV™ PanLassa RDT (Zalgen Labs, LCC, Germantown, MD USA 20876 and Aurora, CO, 80013, Germantown, USA US) – Research for Use Only (RUO)
- Reference standard:** RT-PCR Altona 2.0 kit is used in AE-FUTHA VU laboratory

Setting:

- Ebonyi state: 3 million people, 675,000 people in Abakaliki
- AE-FUTHA (Alex Ekwueme Federal Teaching Hospital, Abakaliki) is a 700 beds tertiary-level hospital
- Supported by MSF since October 2018

Study procedure

- Sample size:
 - Minimum 340 cases, to estimated sensitivity > 90% with 95% certainty
- Patients of **all ages** with suspected Lassa fever + **Informed Consent**
- If yes: RDT will be performed *at bedside*, fingerprick blood
- The RDT was performed by trained health-care staff wearing full Personal Protective Equipment (PPE).
- Visual reading was done twice for each test: at 15 and 25
- Sample for RT PCR taken for routine care – with cycle threshold (Ct value) threshold of 40 to be considered positive
- GPC gene and L gene were the main primers used
- Comparison of RDT from fingerprick and venous sample (RDT repeated in the laboratory)



Source: <https://zalgen.com/product/re-lasv-pan-lassa-antigen-rapid-test-no-50-test-strip-kit/>

Results

- Recruitment during high season 2022-2023
- 217 participants
- Age: median 33 [22.0-44.3]
- Sex: Female: 49.5%; Male 50.5% (**Table 1**)

	PCR positive (N=52)	PCR negative (N=164)*	Total (N=216)*	p value
Sex				
Female	24 (46%)	83 (51%)	107 (50%)	..
Male	28 (54%)	81 (49%)	109 (50%)	0.69†
Age, years	37.5 (22.0-45.0)	32.0 (22.8-43.3)	33.0 (22.0-44.3)	0.70‡

Data are n (%), median (IQR), or p. *One of the 165 participants with negative PCRs had missing data on sex and age. †χ² test with Yates' continuity correction. ‡Wilcoxon rank sum test with continuity correction.

Table 1: Participant characteristics stratified by PCR result

Conclusions


- The Pan-Lassa RDT is **not currently recommended** as a diagnostic or screening tool for suspected Lassa fever cases.
- Marked improvement in **sensitivity and user friendliness is needed** for the RDT to be adopted clinically.
- There remains an urgent need for better Lassa fever diagnostics to promote safety of in-hospital care and better disease outcomes in low-resource settings.

	PCR positive (N=52)	PCR negative (N=165)*
Bedside (cap) RDT at 15 min		
Positive	2 (4%)	0
Negative	47 (90%)	147 (91%)
Invalid	3 (6%)	15 (9%)
Bedside (cap) RDT at 25 min		
Positive	5 (10%)	0 (0%)
Negative	44 (85%)	149 (91%)
Invalid	3 (6%)	14 (9%)
Laboratory (plasma) RDT at 15 min		
Positive	24 (46%)	7 (4%)
Negative	28 (54%)	158 (96%)
Laboratory (plasma) RDT at 25 min		
Positive	26 (50%)	7 (4%)
Negative	26 (50%)	158 (96%)

Data are n (%). *Three participants with a negative PCR result had missing results of the bedside (cap) RDT at 15 min, and two had missing results of the bedside (cap) RDT at 25 min. RDT=rapid diagnostic test.

Table 2: RDT test results stratified by PCR result

- Although the specificity of the Pan-Lassa RDT was high (>90%), sensitivity at bedside using capillary blood was estimated as 4% (95% CI 1–14) at 15 min and 10% (3–22) at 25 min, far below the target of 90%. (**Table 2**)
- The laboratory-based RDT using plasma showed better sensitivity (46% [32–61] at 15 min and 50% [36–64] at 25 min) but did not reach the target sensitivity.
- Among the PCR-positive participants with Lassa fever, positive RDT results were associated with lower cycle threshold values
- Personnel conducting the bedside test procedure reported being hindered by the inconvenient use of full PPE and long waiting procedures before a result could be read.



Source: MSF /

Acknowledgements

We thank all the participants in the RDT study, their relatives, and their caregivers, Abakaliki project staff, both MSF and AE-FUTHA, who contributed to the study. We also thank Anders Boyd, Tony Reid and Armand Sprecher.



	Bedside (cap) RDT at 15 min	Bedside (cap) RDT at 25 min	Laboratory (plasma) RDT at 15 min	Laboratory (plasma) RDT at 25 min
Sensitivity	4.1% (0.5-14.0)	10.2% (3.4-22.2)	46.2% (32.2-60.5)	50.0% (35.8-64.2)
Specificity	100.0% (97.5-100.0)	100.0% (97.6-100.0)	95.8% (91.5-98.3)	95.8% (91.5-98.3)
PPV	100.0% (15.8-100.0)	100.0% (47.8-100.0)	77.4% (58.9-90.4)	78.8% (61.1-91.0)
NPV	75.8% (69.1-81.6)	77.2% (70.6-82.9)	84.9% (79.0-89.8)	85.9% (80.0-90.6)

Data are % (95% CI). Invalid tests were not included. RDT=rapid diagnostic test. PPV=positive predictive value. NPV=negative predictive value.

Table 3: RDT performances of different procedures